# ELKAY Aqua Sentry

Lead, Cyst, Chlorine, Taste& Odor, Atrazine, Lindane, and Particulate Class I Reduction Filter Filter Part No. 56192C For LF2

# FOR COMMERCIAL USE ONLY PERFORMANCE DATA SHEET

NAND CERTIFIES Water Quality Noustry study

LF2 - Aqua Sentry NSF/ANSI 42, NSF/ANSI 53 and NSF/ANSI 372 certified

**IMPORTANT NOTICE:** Read this Performance Data Sheet and compare the capabilities of this unit with your actual water treatment needs. It is recommended that, before purchasing a water treatment unit, you have your water supply tested to determine your actual water treatment needs.

NOTICE: Do not use with water that is microbiologically unsafe or of unknown quality without adequate disinfection before or after the system. Systems certified for cyst reduction may be used on disinfected waters that may contain filterable cysts.

This water filter should be installed and used in compliance with any state or local laws and regulations.

CAUTION: Use on cold water line only. Do not allow installed filter to freeze.

RECOMMENDED USE CONDITIONS					
Capacity	500 Gallons 6 Month Maximum	Maximum Pressure	100 PSI (6.89 bar)		
Flow Rate	0.5 GPM	Minimum Pressure	20 PSI (1 4 bar)		
Temperature	35°-100° F (2°-38° C)		201 01 (111 001)		

NOTE: Performance will vary depending on local water conditions.

### Instructions For Installing Replacement Filters

- 1. Turn off water supply and dispense water to relieve water pressure.
- 2. Mount filter in a suitable location with water flow in direction indicated on filter. This filter is most effective when installed so that water flows through the filter in a downward direction. Do not support the filter using the inlet or outlet tubing. A plastic snap-in bracket, p/n 56190C, is available that can be screwed onto a flat surface.
- 3. Connect the cold water supply to the filter through a shut-off valve and 1/4" O.D. tubing (supplied by the installer). For ease of servicing, the shut-off valveshould be located within easy reach of the filter.
- Note: The inlet and outlet fittings supplied as a part of the filter are John Guest fittings that accept 1/4" O.D. tubing. They may be used with copper, semi-rigid nylon or polypropylene tubing with a rated working pressure equal or greater than encountered in application. Do not use soft braided tubing.
- **4.** Caution: If supply pressure will ever exceed 100 psi, install a pressure regulator to limit the inlet pressure to the filter to 20 100 psi range. To make connections to the filter, push tube end in the fittings past the o-ring to a positive stop approximately 5/8". Ends of tubing must be cut square and free of burrs, sharp ends or scratches that could cut or nick the o-rings or cause a leak.
- 5. Turn water supply on and check fittings for leaks.
- 6. Dispense approximately 1 gallon of water to flush any particulate from the filter.

### LIMITED WARRANTY

The ELKAY Aqua Sentry is warranted to be free from defects in material and workmanship for a period of one year from the date of installation. Warranty is limited to repair or replacement of defective component.

Your Authorized Representative	
rour / autorized representative	
Buyer	

The ELKAY Aqua Sentry filter has been tested according to NSF/ANSI 42 and 53 for reduction of the substances listed below. The concentration of the indicated substance in water entering the system was reduced to a concentration less than or equal to the permissible limit for water leaving the system, as specified in NSF/ANSI 42 and 53.

Substance	Influent Challenge Concentration	Maximum Permissible Product Water Concentration	Reduction Requirements	Minimum Reduction	Average Reduction
NSF / ANSI Standard 42 Aesthetic Effects					
Chlorine Taste & Odor	2.0 mg/L ± 10%		≥ 50%	94.8%	97.4%
Particulates (0.5-<1m)-Class 1	at least 10,000 particles/mL		≥ 85%	99.5%	99.8%
NSF / ANSI Standard 53 Health Effects					
Cysts	Minimum 50,000/L		99.95%	99.95%	99.9%
Lead (pH 6.5)	0.15 mg/L ± 10%	0.010 mg/L		99.3%	99.3%
Lead (pH 8.5)	0.15 mg/L ± 10%	0.010 mg/L		96.6%	96.9%
Atrazine	0.009 mg/L ± 10%	0.003 mg/L		83.2%	91.7%
Lindane	0.002 mg/L ± 10%	0.0002 mg/L		94.5%	98.2%

Testing was performed under standard laboratory conditions, actual performance may vary.

#### **Test Conditions**

Flow Rate	= as noted for filter system		
Inlet Pressure	= 60 psi (4.1 bar)		
рН	7.5 ± 1		
Temperature	$= 68^{\circ}F \pm 5^{\circ}F (20^{\circ}C \pm 2.5^{\circ}C)$		

#### **Operating Requirements**

Pressure	= 20-100 psi (1.4-6.89 bar)	
Temperature	= 35° - 100°F (2° - 37.8°C)	
Turbidity	= 5 NTU Max.	

Recommendation: Filter to be replaced every 6 months or when changes in taste, odor, color and/or flow of the water occurs.

Order Replacement Filter 56192C