



model 7583

Personal Emergency Eyewash (Replacement Bottle)

FEATURES & BENEFITS

IMMEDIATE FLUSHING

The eyewashes provide immediate flushing fluid while the injured person is en route to a primary emergency station.

DUAL EYE PROTECTION

With two ergonomic streams, the DUO eyewash provides simultaneous eye coverage allowing both eyes to be treated with an advanced rinse process for supplementary flushing support.

PH NEUTRAL ADVANTAGE


A sterile pH Neutral phosphate solution is able to neutralize any acidic or alkali substance that enters the eyes to help restore a proper pH balance (only available pH neutralizer on the market). Recommended for chemicals with pH levels above 10 or below 4, though it can also be used for pH neutral contaminants.

ANSI COMPLIANCE

ANSI compliant for supplemental equipment. Reference ANSI Z358.1 Supplemental Equipment/Personal Wash Units (Section 8.1).

OPTIONS

For more information, visit www.hawesco.com or call (888) 640-4297.



SPECIFICATIONS

Model 7583 personal eyewash bottle is a 500 ml (16.9 fl oz) pH Neutral solution which supplies dual eyewash streams to effectively rinse both eyes simultaneously. Recommended wash time is 2 minutes. Fluid carries a 3 year shelf life. Store at room temperature 41-95° F (5-35° C).

APPLICATIONS

Haws Personal Emergency Eyewashes are for use when dirt, dust or other contaminants come into contact with the eyes. Ideal for quick access to first aid in case of eye injury. Mounted at the immediate point of need, this personal eyewash will provide immediate flushing until a primary station can be reached.